

## Product Datasheet

# Recombinant Human Fc gamma RIIA/CD32a, His-Avi (HEK293) (orb1974362)

<b>Description</b>	Recombinant Human Fc gamma RIIA/CD32a, His-Avi (HEK293)
<b>Endotoxins</b>	Not analyzed
<b>Target</b>	FCGR2A
<b>Form/Appearance</b>	Lyophilized or liquid
<b>Concentration</b>	> 1mg/ml
<b>Storage</b>	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
<b>Tag</b>	C-His-Avi
<b>Note</b>	For research use only
<b>Application notes</b>	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
<b>Purity</b>	> 90% as determined by SDS-PAGE
<b>MW</b>	24 kDa
<b>Uniprot ID</b>	<b>P31994</b>
<b>Expiration Date</b>	6 months from date of receipt.