

Product Datasheet

Recombinant SARS-Cov-2 Spike protein(EG.5), (HEK293) (orb1974366)

Description	Recombinant SARS-Cov-2 Spike protein(EG.5), (HEK293)
Endotoxins	Not analyzed
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Tag	C-T4 foldon domain&His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	26.5 kDa
Expiration Date	6 months from date of receipt.