

Product Datasheet

Recombinant SARS-Cov-2 Spike protein(BA.2.86), (HEK293) (orb1974367)

Description	Recombinant SARS-Cov-2 Spike protein(BA.2.86), (HEK293)
Endotoxins	Not analyzed
Target	SARS-CoV-2
Form/Appearance	Liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-T4 foldon domain&10×His
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
MW	136.9 kDa
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.