

Product Datasheet

Recombinant human ALPL protein (orb1974377)

Description	Recombinant human ALPL protein
Target	ALPL
Form/Appearance	Lyophilized or liquid
Concentration	> 1mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	> 90% as determined by SDS-PAGE
Uniprot ID	P05186
Expiration Date	6 months from date of receipt.