

Product Datasheet

Recombinant human PSA protein, C-His (HEK293) (orb1974382)

Description	Recombinant human PSA protein, C-His (HEK293)
Target	KLK3
Form/Appearance	Lyophilized or liquid
Concentration	> 0.5mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	C-His
Note	For research use only
Protein Sequence	18-261/261
Purity	> 90% as determined by SDS-PAGE
MW	28.2 kDa
Uniprot ID	P07288
Expiration Date	6 months from date of receipt.