


Product Datasheet

Human pan CK Ready-To-Use IHC Kit (orb1974483)

Description	orb1974483 is a ready-to-use IHC kit for staining of human pan CK. The kit provides all reagents, from antigen retrieval to cover slip mounting, that require little to no diluting or handling prior to use. Simply apply the reagents to your sample slide according to the protocol and you're steps away from obtaining high-quality IHC data.
Reactivity	Human
Tested Applications	IHC-P
Storage	Please store components at the temperatures indicated on the individual tube labels. The kit is stable for 6 months from the date of receipt.
Note	For research use only
Sample Types	FFPE tissue
Expiration Date	Please enquire.



Immunohistochemical analysis of paraffin embedded human prostate tissue slide using orb1974483 (Human IHC Kit).

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)