

## Product Datasheet

# Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit (orb1974541)

<b>Description</b>	Mouse Glycoprotein Ib Beta Polypeptide, Platelet (GP1Bb) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.2-2000pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Plasma.
<b>Sensitivity</b>	14.7pg/mL
<b>Expiration Date</b>	Please enquire.