

## Product Datasheet

# General Prostaglandin E2 (PGE2) ELISA Kit (orb1974546)

<b>Description</b>	General Prostaglandin E2 (PGE2) ELISA Kit
<b>Reactivity</b>	All
<b>Range</b>	4.68-300pg/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma or other biological fluids.
<b>Sensitivity</b>	1.36pg/mL
<b>Expiration Date</b>	Please enquire.