

## Product Datasheet

# Mouse Diazepam Binding Inhibitor (DBI) ELISA Kit (orb1974550)

<b>Description</b>	Mouse Diazepam Binding Inhibitor (DBI) ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	1.4-90ng/mL
<b>Storage</b>	Store at 4°C & -20°C for 12 months.
<b>Note</b>	For research use only
<b>Sample Types</b>	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
<b>Sensitivity</b>	0.3ng/mL
<b>Expiration Date</b>	Please enquire.