

Product Datasheet

Mouse Beta-Thromboglobulin (bTG) ELISA Kit (orb1974581)

Description	Mouse Beta-Thromboglobulin (bTG) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Storage	Store at 4°C & -20°C for 12 months.
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sensitivity	34pg/mL
Expiration Date	Please enquire.