

Product Datasheet

Goat Anti-Rabbit IgG H&L/PerCP-Cy5.5 (orb1974588)

Description	Goat Anti-Rabbit IgG H&L/PerCP-Cy5.5
Species/Host	Goat
Reactivity	Rabbit
Conjugation	PerCP/Cy5.5
Immunogen	Native Rabbit IgG full length
Preservatives	0.01M PBS (pH7.4)
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.