

## Product Datasheet

# AffiniPure F(ab')<sub>2</sub> Fragment Goat Anti-Rabbit IgG (H+L)/PE (orb1974589)

<b>Description</b>	AffiniPure F(ab') <sub>2</sub> Fragment Goat Anti-Rabbit IgG (H+L)/PE
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rabbit
<b>Conjugation</b>	PE
<b>Immunogen</b>	Native Rabbit IgG
<b>Preservatives</b>	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
<b>Storage</b>	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.