

Product Datasheet

F(ab')₂ Fragment Goat Anti-rabbit IgG H&L / Biotin (orb1974590)

Description	F(ab') ₂ Fragment Goat Anti-rabbit IgG H&L / Biotin
Species/Host	Goat
Reactivity	Rabbit
Conjugation	Biotin
Tested Applications	IHC-P
Immunogen	Native Rabbit IgG
Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	IHC-P=1:400-800
Expiration Date	12 months from date of receipt.