

## Product Datasheet

# SARS-CoV-2 (2019-nCoV) Nucleocapsid Rabbit mAb Antibody (orb1974608)

<b>Description</b>	SARS-CoV-2 (2019-nCoV) Nucleocapsid Rabbit mAb Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, WB
<b>Immunogen</b>	Recombinant SARS-CoV-2 Nucleocapsid Protein (1- 419/419aa)
<b>Preservatives</b>	0.01M PBS(pH7.4)
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>MW</b>	46 kDa
<b>Expiration Date</b>	12 months from date of receipt.