

Product Datasheet

SARS-CoV-2 Antibody (orb1974609)

Description	SARS-CoV-2 Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid Protein (1-419/419aa)
Target	SARS-CoV-2
Preservatives	0.01M PBS(pH7.4)
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	46 kDa
Uniprot ID	P0DTC2
Expiration Date	12 months from date of receipt.