

## Product Datasheet

### CRYAB Antibody (orb1974658)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | CRYAB Antibody  |
| <b>Species/Host</b>        | Rabbit  |
| <b>Reactivity</b>          | Human, Mouse  |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IHC-P   |
| <b>Immunogen</b>           | KLH conjugated synthesised phosphopeptide derived from human alpha B Crystallin around the phosphorylation site of Ser19 FH(p-S)PS      |
| <b>Target</b>              | CRYAB   |
| <b>Preservatives</b>       | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.  |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>MW</b>                  | 20 kDa  |
| <b>Uniprot ID</b>          | <a href="#">P02511</a>  |
| <b>Dilution Range</b>      | IHC-P=1:100-500   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |