

Product Datasheet

phospho-AKT3 (Tyr473)/BF680 Antibody (orb1974989)

Description	phospho-AKT3 (Tyr473)/BF680 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	BF680
Tested Applications	FC
Immunogen	KLH conjugated synthesised phosphopeptide derived from human AKT3 around the phosphorylation site of Tyr473 FS(p-Y)SA
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P31749
Expiration Date	12 months from date of receipt.