

Product Datasheet

CD62p /Biotin Antibody (orb1975202)

Description	CD62p /Biotin Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant mouse P-Selectin protein (42-709/768aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q01102
Expiration Date	12 months from date of receipt.