

Product Datasheet

Phospho-PDPK1 (Tyr9)/BF405 Antibody (orb1975224)

Description	Phospho-PDPK1 (Tyr9)/BF405 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	BF405
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PDK1 around the phosphorylation site of Tyr9 QL(p-Y)DA
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O15530
Expiration Date	12 months from date of receipt.