

Product Datasheet

HBsAg/BF555 Antibody (orb1975306)

Description	HBsAg/BF555 Antibody
Species/Host	Goat
Reactivity	Human, Virus
Conjugation	BF555
Immunogen	Recombinant HBsAg protein
Preservatives	0.01M PBS (pH7.4) with 0.02% Proclin300.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q89714
Expiration Date	12 months from date of receipt.