

Product Datasheet

PACAP receptor/Biotin Antibody (orb1975374)

Description	PACAP receptor/Biotin Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human PACAP receptor (151-250/468aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P41586
Expiration Date	12 months from date of receipt.