

Product Datasheet

phospho-SHC1 (Tyr349)/BF405 Antibody (orb1975446)

Description	phospho-SHC1 (Tyr349)/BF405 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	BF405
Tested Applications	FC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human SHC around the phosphorylation site of Tyr349 HQ(p-Y)YN
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P29353
Expiration Date	12 months from date of receipt.