

Product Datasheet

CXCR3/Biotin Antibody (orb1975538)

Description	CXCR3/Biotin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	WB
Immunogen	Recombinant human CXCR3 protein (1-53/368aa)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O88410
Expiration Date	12 months from date of receipt.