

## Product Datasheet

# phospho-alpha B Crystallin (Ser19)/HRP Antibody (orb1975549)

<b>Description</b>	phospho-alpha B Crystallin (Ser19)/HRP Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse
<b>Conjugation</b>	HRP
<b>Tested Applications</b>	IHC-P
<b>Immunogen</b>	KLH conjugated synthesised phosphopeptide derived from human alpha B Crystallin around the phosphorylation site of Ser19 FH(p-S)PS
<b>Preservatives</b>	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	<a href="#">P02511</a>
<b>Expiration Date</b>	12 months from date of receipt.