

Product Datasheet

ADGRE1/PerCP-Cy5.5 Antibody (orb1975697)

Description	ADGRE1/PerCP-Cy5.5 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	PerCP/Cy5.5
Tested Applications	FC
Immunogen	Recombinant human ADGRE1 (1-530/886aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
Uniprot ID	Q14246
Expiration Date	12 months from date of receipt.