

Product Datasheet

VZV gE/AP Antibody (orb1975702)

| | |
|----------------------------|--|
| Description | VZV gE/AP Antibody |
| Species/Host | Rabbit |
| Reactivity | Virus |
| Conjugation | AP |
| Tested Applications | WB |
| Immunogen | Recombinant Varicella-zoster virus gE protein (31-538/623aa) |
| Preservatives | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | Q9J3M8 |
| Expiration Date | 12 months from date of receipt. |