

Product Datasheet

VZV gE/HRP Antibody (orb1975703)

Description	VZV gE/HRP Antibody
Species/Host	Rabbit
Reactivity	Virus
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant Varicella-zoster virus gE protein (31-538/623aa)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9J3M8
Expiration Date	12 months from date of receipt.