

Product Datasheet

PRCP/HRP Antibody (orb1975741)

Description	PRCP/HRP Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant human PRCP protein (147-496/496aa)
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P42785
Expiration Date	12 months from date of receipt.