

Product Datasheet

ERBB2/HER2/Biotin Antibody (orb1975760)

| | |
|----------------------------|---|
| Description | ERBB2/HER2/Biotin Antibody |
| Species/Host | Rabbit |
| Reactivity | Human |
| Conjugation | Biotin |
| Tested Applications | WB |
| Immunogen | Recombinant human HER2 protein (23-372/1255aa) |
| Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol. |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Uniprot ID | P04626 |
| Expiration Date | 12 months from date of receipt. |