

Product Datasheet

CEBPB Antibody (FITC) (orb1976033)

Description	CEBPB Antibody (FITC)
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CEBP beta around the phosphorylation site of Thr188 + Thr235 PG(p-T)PS
Target	CEBPB
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P17676
Expiration Date	12 months from date of receipt.