

Product Datasheet

BCKDHA, Phospho (S292) Antibody (orb1976077)

Description	BCKDHA, Phospho (S292) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	SW-Size, WB
Immunogen	Surrounding serine 292
Target	BCKDHA, Phospho (S292)
Preservatives	Phosphate Buffered Saline (PBS) with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P12694
Expiration Date	12 months from date of receipt.