

Product Datasheet

TRAP100/MED24 Antibody (orb1976098)

Description	TRAP100/MED24 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Between 939 and C-terminus
Target	TRAP100/MED24
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	O75448
Dilution Range	WB: 1:1000 - 1:5000, IP: 6 µg/mg lysate
Expiration Date	12 months from date of receipt.