

Product Datasheet

XRCC1, Phospho (S518/T519/T523) Antibody (orb1976143)

Description	XRCC1, Phospho (S518/T519/T523) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Surrounding serine 518 and threonine 519 /523
Target	XRCC1, Phospho (S518/T519/T523)
Preservatives	Tris-citrate/phosphate buffer, pH 7 to 8 containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P18887
Dilution Range	WB: 1:1000 - 1:5000, IHC: 1:2000 - 1:10,000. Epitope retrieval with citrate buffer pH 6.0 is recommended for FFPE tissue sections
Expiration Date	12 months from date of receipt.