

Product Datasheet

VSV-G Tag Antibody (FITC) (orb1976152)

Description	VSV-G Tag Antibody (FITC)
Species/Host	Goat
Conjugation	FITC
Tested Applications	ICC, IF
Target	VSV-G Tag
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	Shelf life: 1 year from date of receipt. Storage: 2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	ICC: 1:100 - 1:400, IF: 1:100 - 1:400
Expiration Date	12 months from date of receipt.