

Product Datasheet

MTHFD2L antibody (orb19787)

Description	Goat polyclonal to MTHFD2L
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MTHFD2L
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000. WB: Approx 28kDa band observed in Human Testis lysates (calculated MW of 31.3kDa according to NP_001004346.2). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	KHIQKEIQRGVES
Clonality	Polyclonal
MW	31.3
Entrez	441024
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19787
(0.1ug/ml)
staining of
Human Te...