



Product Datasheet

BMPRI1A Antibody (orb197877)


| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to BMPRI1A |
| Species/Host | Mouse |
| Reactivity | Human, Mouse |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, FC, ICC, IHC, WB |
| Immunogen | Purified recombinant fragment of human BMPRI1A (AA: 179-378) expressed in E. Coli. |
| Preservatives | Ascitic fluid containing 0.03% sodium azide. |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| MW | 60kDa |
| Entrez | 657 |
| Dilution Range | WB: 1/500 - 1/2000; IHC: 1/200 - 1/1000; ICC: 1/200 - 1/1000; RCM: 1/200 - 1/400; ELISA: 1/10000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of PC-3 (Lane 1), ...



Line graph illustrates about the Ag-Ab r...



Flow cytometric analysis of HeLa cells l...