

Product Datasheet

1-epi-Regadenoson (orb1980211)

Description	1-epi-Regadenoson is an α -isomer impurity of Regadenoson, a highly selective adenosine A2A receptor agonist [1].
Storage	-20°C
Note	For research use only
Formula	$C_{15}H_{18}N_8O_5$
Purity	98.00%
MW	390.35
CAS Number	2015222-31-6
SMILES	<chem>CNC(=O)c1cnn(c1)-c1nc(N)c2ncn([C@H]3O[C@H](CO)[C@@H](O)[C@H]3O)c2n1</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)