

Product Datasheet

Leronlimab (orb1980388)

Description	Leronlimab (PRO 140) is a humanized IgG4 anti-CCR5 monoclonal antibody with both anti-HIV and anti-tumor activity, inhibiting CCR-mediated HIV-1 virus and cancer cell metastasis in mouse tumor models. Leronlimab can be used to study HIV, non-alcoholic steatohepatitis (NASH), and breast cancer.
Conjugation	Unconjugated
Target	Leronlimab
Concentration	2.43 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Purity	99.6% (SEC-HPLC)
MW	146.66 kDa
CAS Number	674782-26-4
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)