

Product Datasheet

Acasunlimab (orb1980547)

| | |
|------------------------|--|
| Description | Acasunlimab (GEN1046), a bispecific antibody (bsAb) designed to target PD-L1 and 4-1BB, functions by enhancing both T-cell and NK-cell response via conditional 4-1BB stimulation and continuously inhibiting the PD-1/PD-L1 axis. This compound is applicable in cancer research [1] [2]. |
| Conjugation | Unconjugated |
| Target | Acasunlimab |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Purity | 98.00% |
| CAS Number | 2253937-12-9 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)