

Product Datasheet

Cibisatamab Antibody (orb1980767)

Description	Cibisatamab is a T-cell bispecific antibody that acts as a T-cell splicer. Cibisatamab targets carcinoembryonic antigen (CEA) on cancer cells and CD3 on T-cells, and triggers T-cell killing when the CEA on the surface of the cancer cells ranges from intermediate to high levels. Cibisatamab may be used to study solid cancers such as colorectal cancer.
Conjugation	Unconjugated
Target	Cibisatamab
Concentration	3.46 mg/mL
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Purity	96.5% (SEC-HPLC)
MW	191.1 kDa
CAS Number	1855925-27-7
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)