

## Product Datasheet

### tfap2a antibody (orb19815)

<b>Description</b>	Goat polyclonal antibody to tfap2a
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Zebrafish
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA
<b>Target</b>	Tfap2a (zebrafish)
<b>Preservatives</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Storage</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Note</b>	For research use only
<b>Application notes</b>	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
<b>Protein Sequence</b>	NPNSHSETGSKAGDK
<b>Clonality</b>	Polyclonal
<b>MW</b>	47.1
<b>Dilution Range</b>	ELISA: 1:16000
<b>Expiration Date</b>	12 months from date of receipt.