

Product Datasheet

Secretin (33-59), rat TFA (orb1981592)

Description	Secretin (33-59), rat (TFA), is a 27-amino acid peptide that targets the secretin receptor to enhance the release of bicarbonate, enzymes, and potassium (K+) from the pancreas [1] [2].
Storage	-20°C
Note	For research use only
Formula	$C_{131}H_{217}F_3N_{42}O_{44}$
Purity	98.00%
MW	3141.37
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)