

Product Datasheet

CHRNA1 antibody (orb19820)

Description	Goat polyclonal antibody to CHRNA1
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CHRNA1 / ACHRB
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 55kDa band observed in human brain (cerebellum, hippocampus, frontal cortex) lysates (calculated MW of 56.7kDa according to NP_000738.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QEQEDHDALKED
Clonality	Polyclonal
MW	56.7
Entrez	1140
Dilution Range	ELISA: 1:1000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb19820
(1µg/ml)
staining of
human cere...