

Product Datasheet

MPZ antibody (orb19823)

Description	Goat polyclonal antibody to MPZ
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	Myelin Protein zero
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 28kDa band observed in Rat Spinal Cord lysates (calculated MW of 28.5kDa according to Human NP_000521.1 and of 27.7kDa according to Rat NP_058723.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	DHSRSTKAVSEK
Clonality	Polyclonal
MW	28.5
Entrez	4359
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml
Expiration Date	12 months from date of receipt.



orb19823
(0.1µg/ml)
staining of
Rat Spin...