

Product Datasheet

Ral B protein (orb198260)

Description	Recombinant human Ral B protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IAC
Immunogen	Recombinant GST protein from E. Coli
Preservatives	Provided in 50mM Tris (pH 7.5) solution containing 150 mM NaCl, 10 mM MgCl ₂ , 1 mM DTT, 10% glycerol
Form/Appearance	Provided in 50mM Tris (pH 7.5) solution containing 150 mM NaCl, 10 mM MgCl ₂ , 1 mM DTT, 10% glycerol
Concentration	See vial for concentration
Storage	Product should be stored at -70°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Ral B protein is provided as a functional G-protein and can be used as a Western blot control for antibodies to Ral B protein or in assays to study to protein-protein interactions of Ral B and related proteins with regulatory proteins.
Uniprot ID	P11234
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	6 months from date of receipt.