

Product Datasheet

STING agonist-30 (orb1982600)

Description	STING agonist-30
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{15}H_{16}N_4O_8$
Purity	99.04%
MW	380.31
CAS Number	2951078-67-2
SMILES	<chem>O=C(C1=CC=C([N+])([O-])=O)O1NCCCCNC(C2=CC=C([N+])([O-])=O)O2=O</chem>
Expiration Date	12 months from date of receipt.