

Product Datasheet

JPD447 (orb1983170)

Description	JPD447 is a MAC-0547630 derivative that potentiates UppS inhibitors of β -lactam antibiotics with antimicrobial activity and can be used to study bacterial infections.
Storage	-20°C
Note	For research use only
Formula	$C_{20}H_{23}FN_4$
Purity	99.69%
MW	338.42
CAS Number	2883235-86-5
SMILES	<chem>FC1=CC=C(C=C1)C=2C=NN3C2N=C(C=C3N4CCCCC4)CC</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)