


## Product Datasheet

### BAF312 (Siponimod) (orb198453)

|                        |   |
|------------------------|---|
| <b>Description</b>     | BAF312 is a potent and selective agonist of S1P with EC50 value of 0.39nM for S1P1 receptors and 0.9... |
| <b>Note</b>            | For research use only   |
| <b>Formula</b>         | C <sub>29</sub> H <sub>35</sub> F <sub>3</sub> N <sub>2</sub> O <sub>3</sub>                            |
| <b>MW</b>              | 516.59501   |
| <b>CAS Number</b>      | [1230487-00-9]  |
| <b>SMILES</b>          | N(CC1=CC=C(/C(=N/OCC2=CC=C(C3CCCCC3)C(C(F)(F)F)=C2)/C)C=C1CC)1CC(C(O)=O)C1                              |
| <b>Expiration Date</b> | 12 months from date of receipt.   |



Chemical structure of BAF312 (Siponimod)