

Product Datasheet

cIAP1 Ligand-Linker Conjugates 16 (orb1984635)

Description	cIAP1 Ligand-Linker Conjugates 16 refers to a specific type of ligase ligand-linker conjugate, possessing the ability to engage with E3 ligases, which finds application in the synthesis of PROTACs.
Storage	-20°C
Note	For research use only
Formula	$C_{29}H_{37}N_9O_5S$
Purity	98.00%
MW	623.73
CAS Number	2143956-63-0
SMILES	<chem>Cc1[nH]nc(Nc2ncnc3cc(OCCCN4CCN(CC4)c4cnc(cn4)C(O)=O)c(cc23)S(=O)(=O)C(C)(C)C)c1C</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)