

Product Datasheet

Lasmiditan (COL-144; LY573144) (orb198469)

Description	Lasmiditan (COL-144; LY573144) is a high-affinity, highly selective 5-HT1F receptor agonist(Ki=2.1 n...
Note	For research use only
Formula	$C_{19}H_{18}F_3N_3O_2$
MW	377.3672
CAS Number	[439239-90-4]
SMILES	<chem>FC1=CC(F)=C(C(NC2=CC=CC(C(C3CCN(C)CC3)=O)=N2)=O)C(F)=C1</chem>
Expiration Date	12 months from date of receipt.